

# microtech™ QuickTIPS

## Steps for Success—An Initial Fitting Protocol


### Product Availability:

Propel, Propel Tinnitus.

### Required Inspire

1. Connect hearing aids to the fitting computer via the appropriate programming accessories
2. **Auto Path** launches automatically when hearing aids that have not been previously read are detected by Inspire
3. Ensure all icons are **checked** and highlighted in **orange**
4. Click **Next** to start the routine
5. Follow all onscreen prompts
6. Make sure Receiver Matrix and Acoustic Options selection in Inspire must match device configuration to accurately reflect the frequency response
7. Choose the **Voice Language** for Indicators
8. Click **Next**
9. Select **Experience Level**
10. Click **Close** to exit Auto Path and open Quick Fit


- Programming cables and connectors spreadsheet is available via Inspire Online
- Auto Path protocol varies based on product family and technology level
- Select Experience Level 3 for patients experienced with amplification
- Select Experience Level 1 for patients new to amplification



Inspire Auto Path

### QuickFit


11. Use the controls in the bottom center of the QuickFit screen to adjust for patient comments regarding output, gain and occlusion
12. Click on the **User Control Configuration** icon to enable the user control as memory, volume control, or both. (Propel Tinnitus also offers an adjustment for the tinnitus level)
13. Enable **Memory 4** as the preferred automatic telephone response via the memory tabs in lower corner of the frequency response screen



QuickFit

### Speech Mapping

14. Click on **Speech Mapping** in the left navigation bar
15. Place a check for soft (50) and loud (80) Target and Response curves in the Display menu
16. Select the movie camera icon to capture a 10 second video using live or recorded voice



Speech Mapping


# microtech™ QuickTIPS

17. Ensure the top capture is between the soft and loud curves to demonstrate audibility and comfort

## FrequencyTrax


18. Click on **FrequencyTrax** on the left navigation bar
19. **Bandwidth** adjusts the Source and Target regions to maximize speech clarity
20. **Gain** adjusts the amplification of the recreated cues to maximize audibility (ten 1dB increment steps)
21. **Expert Assistant** accessed via the left navigation bar or QuickFit screen included additional fine tuning assistance
22. Click on the **Memory Tabs** in the lower corner of the frequency response graphs to enable or disable FrequencyTrax per memory

FrequencyTrax in Inspire



## Indicators


23. Select **Indicators** from the left navigation bar
24. Click on an Indicator icon to highlight it to verify audibility
25. Use the **Play** icon to test audibility through the hearing aid
26. Use the **Speaker** icon to demo the indicator through the computer speakers
27. Use the **+** icons to increase or decrease presentation level as needed
28. **Languages** may be adjusted from this screen
29. Available **Voice Languages** may be demonstrated via the Media File Player



Indicators Screen

## Fitting Summary

30. Select Fitting Summary from the left navigation bar
31. Hover over any screen label to highlight the hyperlink
32. Click on the hyperlink to navigate to the adjustment screen if needed
33. Enter detailed Fitting comments and save or print to the patient's file if required.



Fitting Summary